

## **BLOOD DNA EXTRACTION – MINI KIT**

### **Kit Components**

<b>Components</b>	<b>Storage</b>
Blood Lysis Buffer 1	RT
Solution A	RT
Blood Lysis Buffer 2	RT
Proteinase K	2 - 8°C
Proteinase K Buffer	2 - 8°C
Blood MagNa Mix	RT
Blood Wash Buffer 1	RT
Blood Wash Buffer 2	RT
Blood Elution Buffer	2 - 8°C
MagNa Stand (optional)	RT

### **Materials needed prior to start**

- Add 96 - 100% Ethanol to Wash Buffers as indicated on the bottle.
- Set the water bath/heat block at 56°C.
- Reconstitute Proteinase K with Proteinase K Buffer and store at 2 - 8°C.

## **PROTOCOL**

### **Sample**

- Human blood stored in EDTA\Heparin\Citrate\Fluoride treated vacutainer (Fresh, old or frozen)

### **Blood Lysate Preparation**

1. To a sterile 1.5ml microcentrifuge tube, add **400 µl of blood sample**.
2. Add **1 ml of Blood Lysis Buffer 1** and mix the contents by inverting the tube 4 - 5 times.
3. Add **50 µl of Solution A** and mix the contents by inverting the tube 4 - 5 times (*do not pipette mix*).
4. Incubate at room temperature for 5 minutes.



5. Centrifuge at **8000 rpm for 2 minutes at room temperature.**
6. Carefully discard the supernatant and add **1 ml of Blood Lysis Buffer 1.**
7. Resuspend the pellet in the buffer (*thoroughly by pipette mixing*).
8. Centrifuge the sample at **8000rpm for 2 minutes** at room temperature and carefully discard the supernatant.
9. To the pellet, add **500 µl of Blood Lysis Buffer 2** (*completely resuspend the pellet by pipette mixing*).
10. Add **20 µl of Proteinase K** and mix the contents by pipetting.
11. Incubate at **56°C for 10 minutes.**
12. Transfer the lysate completely to a fresh 1.5 ml microcentrifuge tube.

#### **DNA Binding & Washing**

*Note: Vortex the Blood MagNa Mix to ensure complete dispersion of the particles.*

13. To the lysate, add **350 µl of Blood MagNa Mix.** Gently mix the contents by inverting the tube for 10 - 12 times. Incubate at room temperature for 5 minutes.
14. Place the tube on MagNa Stand for 5 minutes.
15. Carefully discard the supernatant without removing the tube from MagNa Stand (*make sure that the pellet is not disturbed*).
16. Add **500 µl of Blood Wash Buffer 1,** remove the tube from MagNa Stand and resuspend the pellet by pipette mixing for about 5 - 8 times (*ensure complete dispersion of the particles*).
17. Place the tube back on MagNa Stand until the solution becomes clear (*30secs - 1min.*).
18. Carefully discard the supernatant without removing the tube from MagNa Stand (*make sure that the pellet is not disturbed*).
19. Add **500 µl of Blood Wash Buffer 2** and gently invert mix the tube without removing from MagNa Stand for 5 - 6 times to wash the pellet (*Surface wash only*).
20. Discard the supernatant without removing from the MagNa Stand.
21. Repeat step 19-20.



22. Air dry the pellet without removing the tube from MagNa Stand at room temperature for 10 - 15 minutes.

*Note: Do not over dry the pellet.*

### **DNA Elution**

23. After drying, remove the tube from MagNa Stand and add **100 µl of Elution Buffer**.
24. Resuspend the pellet by pipette mixing 10 - 12 times (*ensure complete dispersion of particles*).
25. Incubate at **56°C for 5 minutes** with intermittent tapping.
26. Place the tube on MagNa Stand for 5 minutes or until the solution gets cleared.
27. Carefully transfer the extracted DNA to a sterile microcentrifuge tube, without removing the tube from MagNa Stand. Make sure that the pellet is not disturbed.
28. Discard the MagNa particles.

*Note: In the elution step, if the MagNa particles take more than 10 minutes for clearing, then spin the tubes at 14,000 rpm for 5 minutes and place the tubes back on MagNa Stand until solution becomes clear. Transfer the supernatant to a fresh tube.*

### Downstream applications

DNA obtained from this method is ready to use for any downstream application without any further precipitation or purification step.

### Troubleshooting Guide

Observation	Possible Causes	Suggested Solution
<b>Low DNA yield or Poor Quality</b>	Incomplete Lysis	Use suggested amount of Proteinase K for the specified time.
		Ensure proper re-suspension of pellet while adding the lysis buffers.
		Make sure that the Lysis incubation temperature and time is followed correctly.
		In some cases when DNA content is very high, the lysate may appear slimy. Ensure thorough pipette mixing during and after lysis, followed by the addition of MagNa Mix in such circumstances.
	Incorrect reagent volumes were used	Use the exact volumes of reagents mentioned in the protocol.
	MagNa Mix was improperly handled	Resuspend the MagNa Mix by vortexing prior to use.
	MagNa Mix was disturbed or lost during binding or washing steps	Carefully remove the supernatant from the tube without removing the tube from the MagNa Stand without disturbing the MagNa particles.
Improper elution	Completely resuspend the MagNa particles in elution buffer before incubation at 56°C	

		for elution.
	Ethanol is not added to wash buffers	Add 100% ethanol to wash buffers prior to use as mentioned on the bottles.
<b>Poor performance of extracted DNA in downstream applications</b>	Ethanol carryover	Air dry the MagNa particles properly after washing steps, to remove the ethanol completely, but do not over dry the pellet.
	Salt carryover	Ensure that the correct amount of ethanol is added to the Wash Buffers and two wash steps are performed with Wash Buffer 2.

*\*This kit is designed to extract Genomic DNA from blood clot and PBMCs as well. For protocols & detailed information, please visit [www.maggenome.com](http://www.maggenome.com)*